


Aloe Vera Gloves

Healthy & Soft Hands

Aloe gloves are clinically proven to give biologists softer, healthier and cleaner hands.



Onset of skin quality improvement with Aloe Gloves

West & Zhu, 2003
Am. J. Infect. Control.
2003; 31: 40-2



Aloe exam gloves were developed specifically for those of you who change gloves and wash your hands frequently throughout the day. Aloe exam gloves are coated on the inside with aloe vera gel, working to soothe and hydrate your skin every time you wear them. The aloe vera used in the gloves is pure and pasteurized, and contains no oils or grease.



D-MARK BIOSCIENCES
E-mail: info@d-markbio.com
Toll Free: 1-800-784-7485